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| TEST ITEM                          | UNIT  | SAMPLE | RESULT       | TEST METHOD   |
|------------------------------------|-------|--------|--------------|---------------|
| Residual solvent : Divinyl sulfone | mg/kg | -      | Not Detected | By the client |

\* GC/MS : Gas chromatograph mass spectrometer  
 - Test period : Mar. 16  
 - Assesment : Package of ...  
 - Usage of Residual Solvent

NOTE : 1. The test result is for the client and do not use for other purpose.  
 2. This test report shall be used only within the purpose of its defined usage and shall not be used for public relation, advertisement and lawsuit.  
 3. This test report is only valid when printed on KTR original report paper with hologram and when re-issued by KTR. The copy and the electronic file of the test report are only for reference.

| TEST ITEM                          | UNIT  | SAMPLE | RESULT       | TEST METHOD   |
|------------------------------------|-------|--------|--------------|---------------|
| Residual solvent : Divinyl sulfone | mg/kg | -      | Not Detected | By the client |

2017.03.16.  
**Korea Testing & Research Institute**  
 President *Byun, Jong-Ryeop*

1 of Total : 1 Page(s)

KTR-QP-109-F01-06(03) AIZB 8.78(1)

**KTR** KOREA TESTING & RESEARCH INSTITUTE

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## Odniesienia

- Stern R. Hyaluronan catabolism: a new metabolic pathway. *Eur J Cell Biol.* 2004;83:317-25. [PubMed]
- L. Ambrosio, A. Borzacchiello, P. A. Netti, and L. Nicolais. "Properties of new materials: rheological study on hyaluronic acid and its derivative solutions." *Journal of Macromolecular Science-Pure and Applied Chemistry.* vol. 36. no. 7-8, pp. 991-1000, 1999. [View at Google Scholar](#) [View at Scopus](#)
- G. D. Monheit and K. M. Coleman. "Hyaluronic acid fillers." *Dermatologic Therapy.* vol. 19. no. 3. pp. 141-150. 2006. [View at Publisher](#) [View at Google Scholar](#) [View at Scopus](#)
- R. D. Price, M. G. Berry, and H. A. Navsaria. "Hyaluronic acid: the scientific and clinical evidence." *Journal of Plastic, Reconstructive & Aesthetic Surgery.* vol. 60, no. 10. pp. 1110-1119, 2007. [View at Publisher](#) [View at Google Scholar](#) [View at Scopus](#)
- S. L. Matarasso, J. D. Carruthers, and M. L. Jewell. "Consensus recommendations for soft-tissue augmentation with nonanimal stabilized hyaluronic acid (Restylane)." *Plastic & Reconstructive Surgery* vol. 117, no. 3, pp. 3S-34S, 2006. [View at Publisher](#) [View at Google Scholar](#) [View at Scopus](#)
- A. Tezel and G. H. Fredrickson. "The science of hyaluronic acid dermal fillers." *Journal of Cosmetic and Laser Therapy,* vol. 10. no. 1, pp. 35-42, 2008. [View at Publisher](#) [View at Google Scholar](#) [View at Scopus](#)
- U. B. G. Laurent and R. K. Reed. "Turnover of hyaluronan in the tissues," *Advanced Drug Delivery Reviews,* vol. 7, no. 2, pp. 237-256, 1991. [View at Publisher](#) [View at Google Scholar](#) [View at Scopus](#)
- A. Almond, "Hyaluronan," *Cellular and Molecular Life Sciences.* vol. 64, no. 13. pp. 1591-1596, 2007. [View at Publisher](#) [View at Google Scholar](#) [View at Scopus](#)
- L. Mayol, M. Biondi, L. Russo, B. M. Malle, K. Schwach-Abdellaoui, and A. Borzacchiello. "Amphiphilic hyaluronic acid derivatives toward the design of micelles for the sustained delivery of hydrophobic drugs," *Carbohydrate Polymers.* vol. 102, no. 1, pp. 110-116. 2014. [View at Publisher](#) [View at Google Scholar](#) [View at Scopus](#)

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Company name : BioPlus Co., Ltd.